

SOFT

Wrist blood pressure monitor



SOFT INFLATE TECHNOLOGY:
More comfortable, more precise



Average of last 3 measurements



Irregular heart rhythm indicator



Hypertension level indicator



5 years warranty

The new Medel SOFT wrist blood pressure monitor comes with an innovative measuring technology able to detect the measurement during the inflation, for more comfortable and accurate measurements. A scale gives you the pressure values indication (position 1-3 means an eventual hypertension). Medel SOFT also comes with a very large display, so you can read the values in a more visible way, also checking if an eventual irregular heart rhythm has been detected during the measurement. You can also check the average of the last 3 values. 5 Year Warranty. Ila Class Medical Device.



Movement sensor



Wrist circumference 13.5-21.5




Battery level indicator



60 memory positions



Batteries included



TECHNICAL SPECIFICATIONS

Technology	Oscillometric
Accuracy	pressione + / - 3 mmHg battiti + / - 5%
Measurement range	SIST: da 60mmHg a 230 mm Hg DIAST: da 40mmHg a 130mmHg PULS: 40-199 battiti/min
Power supply	2 Alkaline AAA
Memory	60
Users	1
Date and time	✓
Motion sensor	✓
LCD size (bxh)	3,2 cm x 4,5 cm
Cuff tipe	wrist - soft technology
Cuff size	Universal 13,5 - 21,5 cm



ACCESSORIES

- Storage case, universal cuff

95216

